

## HIV 101: FDA-Approved Rapid HIV Antibody Screening Tests

Assay by Specimen Type*	Manufacturer	Sensitivity (95% CI)	Specificity (95% CI)	Window Period for Result Validity†
<b>Whole Blood (finger stick or venipuncture)</b>				
Clearview HIV 1/2 STAT-PAK	Inverness Medical Professional Diagnostics (www.invernessmedicalpd.com)	99.7% (98.9-100)	99.9% (99.6-100)	15-20 minutes
Clearview COMPLETE HIV 1/2	Inverness Medical Professional Diagnostics (www.invernessmedicalpd.com)	99.7% (98.9-100)	99.9% (99.6-100)	15-20 minutes
OraQuick ADVANCE Rapid HIV-1/2 Antibody Test§	OraSure Technologies, Inc. (www.orasure.com)	99.6% (98.5-99.9)	100% (99.7-100)	20-40 minutes
Uni-Gold Recombigen HIV Test	Trinity Biotech (www.unigoldhiv.com)	100% (99.5-100)	99.7% (99.0-100)	10-12 minutes
<b>Serum or Plasma</b>				
MultiSpot HIV-1/HIV-2 Rapid Test	BioRad Laboratories (www.biorad.com)	100% (99.9-100)	99.9% (99.8-100)	Immediately to up to 24 hrs
Reveal G3 Rapid HIV-1 Antibody Test	MedMira, Inc. (www.medmira.com)	99.8% (99.0-100)	98.6% (98.4-98.8)	Must be read immediately

\* When tests may use other specimen types, it is listed as a table footnote.

† As measured from last step of testing process

\_ Can also be used with a serum or plasma sample

§ Can also be used with an oral fluid specimen or plasma sample

Data adapted from References 11 and 12 of the Main Article.